



MICROSCOPE SLIDE WARMER

Code no: 06/0045

WARMS MICROSCOPE SLIDES FOR SEMEN MOTILITY EXAMINATION

An economical alternative to a heated stage, the Rotech Slide Warmer will warm a handful of slides to 35°C for motility examination of semen.

Semen motility examination must be carried out using a warm slide, because low temperature will have a negative effect on motility, leading to the possibility of an inaccurate motility score being given.

A heated stage will keep a slide warm during the examination, but it is not practical for warming it prior to preparing the sample.

The **Rotech Slide Warmer** simply needs to be plugged into the mains supply and then slides placed onto its surface. The thermostatically controlled device will then **maintain a temperature of 35°C**, so that when you come to place a drop of semen for examination, there will be no temperature shock.

And because motility checking is so quick, you may not even need to have a heated stage.

Specification:

Dimensions: 155 x 155 x 25 mm

Power Requirements: 220v or 110v AC (*please specify when ordering*)

Surface Temperature: Pre-adjusted to 35°C.



Rotech Livestock Equipment Limited

Woodhorn Lane, Oving, Chichester
PO20 2BX, England

Tel: (+44) (0)1243 787115

Fax: (+44) (0)1243 531847

E-mail: admin@rotechuk.com

Website: www.rotechuk.com

